

Product Name : Drag Co Efficient Apparatus	Product Code : NLAB-TECHNICALAB32002
Description :	
Drag Co Efficient Apparatus	
Technical Specification :	
<p>The apparatus has been designed to introduce students to the fundamental characteristics of the behavior of the particle system, in particulars the relationship between the drag coefficient of falling particles and their Reynolds number value. Particles covering a range of sizes and densities are supplied and the experiments are conducted by allowing single particles to fall through a number of different liquids contained in vertical glass tube. The rate of fall of the particle is determined by timing their passage between two marks on the walls of the glass tubes. Observation of particle movement is aided by the provision of a shielded fluorescent tube light mounted on the backboard. Particles can be removed from the bottom of the tubes without the necessity of draining the liquid. A valve system is provided at the bottom of each tube to allow the particles to be removed with the minimum loss of liquid.</p> <p>Experiment:</p> <ul style="list-style-type: none">To measure the drag co-efficientTo plot the graph between drag co-efficient vs. Modified Reynolds numberTo study the motion of a solid particle through a liquid <p>Utilities Required:</p> <ul style="list-style-type: none">Electricity supply 100 watts, 220 , 1 phaseFloor Area 0.75 m x 1.5 mLiquids of different viscosity	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003