

<b>Product Name :</b> Extraction Thimbles	<b>Product Code :</b> NLAB-CHEMICALSIND19007
<b>Description :</b>	
Extraction Thimbles	
<b>Technical Specification :</b>	
<p>Thimbles are used in soxlet apparatus for dust sampling, and for extracting soluble components using organic solvents. These are typically employed, where multiple extraction cycles are required to dissolve the desired compound. Extraction Thimbles are often used in Soxhlet extractors or automated extractors. Thimbles are made of cellulose, glass fiber, or silica fiber. Typical applications include analysis of fats, oils, grease, organics, pesticides, pollutants, and additives in plastics/rubber materials. Extraction units equipped with extraction Thimbles minimize the hands-on time and ensure a complete extraction.</p>	

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003