

Product Name : Haemacytometer Glass	Product Code : NLAB-EDUSCHOOLEQ173034
Description :	
Haemacytometer Glass	
Technical Specification :	
<p>Each plateau features enhanced Neubauer rulings, the ruled surface being 0.1 mm below the cover glass, limiting the volume of blood or fluid over a square mm to 0.1cu.mm and over each of 400 squares (within the central square mm) to 0.00025 cu. mm.</p> <p>These Bright Line glass Haemacytometers are designed for accurate counting of cells in blood or other fluids, with one-piece glass construction, an H-shaped moat to form two counting areas (plateaus) and a "V" slash at the loading side of each plateau.</p> <p>Contact of the flat, polished cover glass surfaces with cover glass supports produces an exact volume of fluid over the counting area.</p> <p>Supplied as a double chamber with two cover glasses or a kit which includes a double counting chamber, cover glasses and red and white pipettes, all in a protective case.</p> <p>The difference in surface tension characteristics between the metallic surface on the chamber and the polished cover glass assures smooth capillarity for precise loading and more even cell distribution.</p> <p>Spare/replacement glass cover glasses are also available.</p> <p>Improved Neubauer ruling pattern</p> <p>Bright Line (Rhodium coated) counting chamber</p> <p>Volume: 0.1 µL</p> <p>Cell Depth: 0.100 mm ±2 %</p> <p>To clean the counting chamber after use: remove the cover glass and clean with water or a mild cleaning solution (e.g. 10% bleach solution). Dry the counting chamber with a soft cloth/wipe, or rinse with acetone.</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003