Phone: +91-9896600003

Email: sales@naugralabequipments.com

Product Name:

ICU Ventilator

Product Code:

NLAB-SAFETYLABP11022



Description:

ICU Ventilator

Technical Specification:

An Optimal Combination of Invasive and Noninvasive Ventilator.

As noninvasive ventilation, it is used increasingly in a wide range of clinical situations, we offer a dual solution AN359 which combines the advantages of a flexible noninvasive ventilator with a full-featured invasive ventilator for the ICU

Ventilation Modes:

 $\label{eq:cv} VCV(A/C), PCV(A/C), PRVC, SIMV(VCV) + PSV, SIMV(PCV) + PSV, SIMV(PRVC) + PSV, SPONT/CPAP + PSV, BIVENT + PSV, NIV/CPAP, NIV-T, NIV-S/T$

T0da1006 1mme:	
Respiration Rate:	
Tiasps	
Tslope:	
Մթակաse:	
It: El (Ratio):	

Phone: +91-9896600003

Email: sales@naugralabequipments.com

Pi@⊈sure (-20~0 cmH2O, above PEEP), **Flow** (0.5~20 LPM)

TriggemSH20sitivity:

PEEP:

Psupportt20

PinspcmH2O

Alarm:

Paw high/low MVe high/low Circuit disconnnect

FiO2 high/low Inspiration/Expiratory tidal volume low

High Respiration Rate Apnea AC Failure Nebulizer On

Low Battery Air /O2 supply down High/Low PEEP

Leakage out of range Occlusion

Siza Entre Color touch screen (detachable)

5040 p01 y2 65-a0s 6 MPa

PANDIANOS 24400 Ny; 50 Hz/60 Hz

Com 200 20 Infortation in the contract of the

BRADENSING (WXDXH) 366 mm (Main Unit), 547 mm x

675 mm x 950 mm (Cart)

Weight (Main Unit), 25 kg (Cart)

Superior Mobile ICU ventilator

Comprehensive ICU ventilator including BIVENT and PRVC

Compact, big capacity battery, no air compressor, intra-hospital mobility

Flexible device configuration: equipped on a trolley, bed or ceiling pendant

Cost Effective Solution:

Unique metal-based, autoclavable, heated exhalation valve

Built-in flow sensor, non-consumable design

Upgradeable ventilation system software, with an available USB port

Optimal patient-ventilator synchrony, increase patient comfort

The Unique Leak Compensation System - Keep precise control on the tidal volume of each breath delivered to the patient by adjusting compensation dosage automatically

Phone: +91-9896600003

Email: sales@naugralabequipments.com

Advanced Trigger Technique - Enhance sensitivity, avoid spurious triggering

Safe Ventilation Through Whole Treatment Phase Initial Treatment Phase

Noninvasive ventilation mode associated with decreased intubation rates, shortened patient stays, improved patient comfort, and a reduced risk of cross infection

Preset patient's height and IBW. Reduce clinician's workload

Stable Condition Phase

PRVC and BIVENT employ lung-protective strategies, delivering intelligent ventilation

Comprehensive lung mechanics monitoring include compliance, airway resistance, PEEPi and time constant

Three waveforms & three loops with user-friendly display provide a continuous monitoring of the patient's condition

Weaning Phase

Various ventilation modes enhance the weaning process

The unique trigger and leakage compensation system safeguards each and every patient breath resulting in smooth and comfortable breathing, avoiding extra workload on the patient and promoting recovery

RSBI and WOB provide accurate reference for weaning

Rehab Phase

Data export port provides connection to hospital monitors and Patient Data Management Systems

Provides pressure support for the patient when spontaneous breathing is present

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com **Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003