Phone: +91-9896600003

Email: sales@naugralabequipments.com

## **Product Name:**

Naugra Lab Flocculator Jar Test Apparatus

## **Product Code:**

NLAB-LABORATORYEM64001



## **Description:**

Naugra Lab Flocculator Jar Test Apparatus

## **Technical Specification:**

Flocculator Jar Test Apparatus

## Features:-

Accepts two, four & six samples

Digital indication of speed and time

Two pre-set programs

Phone: +91-9896600003

Email: sales@naugralabequipments.com

uilt in illumination for easy viewing of samples			
Flocculation/Jar Test Apparatus, is a flocculation (jar) tests on water and effluent samples.			
A digital display clearly indicates the speed of the rotational stainless steel paddles, which can be varied from 25 to 250rpm.			
The digital timer can be set to count down from 1 to 99 minutes.			
After count down, the paddles stop and an alert sounds.			
The timer then counts up in minutes in order to measure settling time.			
Two pre-set programs allow routine speed and time parameters to be input and then selected by simply pressing one button.			
Parameters are stored in the memory and can be recalled or altered at any time.			
The intuitive touch sensitive keypad is totally flat for easy cleaning.			
Specifications:-			
<b>B</b> lo. of Jar			
25 de 200 nge			
Pre-set programs			

Phone: +91-9896600003

Email: sales@naugralabequipments.com

Dtal Display	<b>īito⊕</b> 9minutes	
	Counctro Issel Ke ij tipvæd	
Weir50Hz, 100W	Dīgūtal Display	
	<b>?‰</b> k;r50Hz, 100W	

# Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com **Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003