

**Product Name :**  
Particle Size Analysis Set Pipette Method

**Product Code :**  
NLAB-SOILELAB140041



**Description :**

Particle Size Analysis Set Pipette Method

**Technical Specification :**

This is for the determination of the sub sieve particle distribution in a soil sample by mechanical analysis. An analysis of this kind expresses quantity the proportions by weight of the various sizes of particles present in the soil. It is recommended as a standard procedure to use dispersion agent to avoid flocculation. The apparatus consists of a sliding panel which moves up and down by means of a screw allowing Anderson pipette fixed to it to be raised or lowered vertically. A sedimentation tube is held by a laboratory clamp provided on the stand below the pipette. The depth of immersion is measured by a scale graduated in mm at the side of the sliding panel. Supplied complete with Anderson pipette 10ml. At the side capacity made from glass, and a sedimentation tube also of glass of 500ml capacity and 50nos. Test form pads.

Particle size analysis set: Pipette method

Accessories & Spares: (on extra cost)

Sedimentation Pipette (Anderson pipette) 25ml.

Sedimentation tube 100ml

Sedimentation pipette 10ml

Sedimentation tube 500ml

Test forms pad of 50

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003