Phone: +91-9896600003

Email: sales@naugralabequipments.com

Product Name : RTD Studies in Packed Bed	Product Code : NLAB-EDUCATIONTVET16001
Description :	
RTD Studies in Packed Bed	
Technical Specification :	
RTD Studies in Packed Bed -	
Liquid is fed to the reactor through liquid distributor, fitted at the bottom of the	ne column.
Rota meter is provided to measure the flow.	
Treata meter le previdea le mededie die mem	
The setup consists of a glass column packed with Rasching Rings.	
For understanding the R.T.D. characteristics, a special arrangement to	
inject tracer into the lower end of reactor (for liquid and gas both) is provided.	
F	
The flow rate can be adjusted by operating the needle valve provided on Ro	ota meter.
Samples can be taken periodically from the top outlet of reactor Experiments	s.

Phone: +91-9896600003

Email: sales@naugralabequipments.com

To plot RTD curve for Plug Flow Reactor.

To determine the Dispersion No Utilities Required Water Supply & Drain Instruments, Laboratory Glassware and Chemicals required for analysis as per the system adopted.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com **Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003