Phone: +91-9896600003

Email: sales@naugralabequipments.com

| Product Name : Thermal Conductivity of Metal Rod | Product Code : NLAB-EDUCATIONTVET16017 |
|---|---|
| | |
| | |
| | |
| | |
| | |
| | |
| Description : | |
| Thermal Conductivity of Metal Rod | |
| Technical Specification : | |
| Thermal Conductivity of Metal Rod - | |
| The experimental set up consists of a metal bar, one end of which is heated by an electric heater while the other end of the bar projects inside the cooling water jacket. By varying the heat input rates, data can be obtained. Heat Input to the heater is given through variac. | |
| The temperature of the bar is measured in different sections. | |
| A cylindrical shell filled with insulating material surrounds the middle portion | n of the bar. |
| Water at a constant rate is circulated through the jacket and its flow rate an | nd temperature rise are noted. |
| | |
| | |
| | |

Phone: +91-9896600003

Email: sales@naugralabequipments.com

| Experimentation:- |
|---|
| |
| |
| |
| To plot the temperature distribution along the length of Bar. |
| |
| |
| To determine the thermal conductivity of given bar at various temperatures. |
| To determine the thermal conductivity of given bar at various temperatures. |

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com **Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003