

Product Name :
V.D.R.L. Rotator

Product Code :
NLAB-SCIENTIFICEQ36005



Description :

V.D.R.L. Rotator

Technical Specification :

APPLICATIONS OF V.D.R.L. ROTATOR

Venereal Disease Research Laboratory (V.D.R.L) Rotator is used for V.D.R.L. tests, blood-grouping tests, and various agglutination tests or for mixing of solutions in small bottles, flasks or beakers.

Useful for serological Test and other tuned clinical diagnostic procedures requiring gentle agitation.

SALIENT FEATURES & CONSTRUCTION OF V.D.R.L. ROTATOR

It is a compact bench top model unit

The shaker is mounted on ball bearing crankshafts to give an orbital circular motion to the flasks.

The rotator has covered metal platform Size 130 x 130 with spring, holder clamp.

The diameter of the orbit is approximately – 40 mm

Compact counter balanced drive mechanism ensures high stability and reliability even in continuous operation & uneven load distribution.

Speed range is 70 – 210 RPM.

High Speed (Model) range is up to 350 RPM.

Fitted with heavy duty motor having variable speed.

Speed can be controlled through variable speed regulator knob.

Control panel plate fitted on the front having indicator lamp ON/OFF switch & speed regulator knob.

Supplied complete with cord & plug to work on Single Phase 220 / 230 volts A.C. supply.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003