

Product Name :
250/25 Kn Automatic Compression Press

Product Code :
NLAB-EDUCATIONTVET15052



Description :

250/25 Kn Automatic Compression Press

Technical Specification :

Automatic charge cycle

Display and control by integrated 7 "color touch screen, protection class IP31

Memory for 100 test profiles and 6 favorites. 2Gb of memory allows the storage of more than a million tests

LAN interface

Comes with download software for Windows

Accuracy better than +/- 1% in the calibrated range

Two test spaces

Frame 250 kN:

Horizontal passage max 225 mm

150 mm upper and lower shelves

Frame 25 kN:

Vertical passage max 205 mm

150 mm upper and lower shelves

Piston stroke 15 mm

Supplied complete with accessories for compression and bending of 4x4x16 cm specimens

Set of spacers

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003