

Product Name :
Adiabatic Batch Reactor

Product Code :
NLAB-EDUCATIONTVET16031

Description :

Adiabatic Batch Reactor

Technical Specification :

Adiabatic Batch Reactor -

This reactor is operated under conditions like adiabatic (Temperature of Reaction Mass keep on changing), Perfectly mixed (composition of the reaction mixture is uniform throughout), and constant volume (volume of the reaction mixture within the reactor remains constant, there is no appreciable change in the density of reaction mass). An adiabatic batch reactor is a closed system with no input and output streams. This set-up is used to study a catalytic homogeneous reaction under adiabatic condition.

The set-up consists of a Dewar flask of capacity 1 lit. fitted with a stirrer.

The flask is double insulated by ceramic wool insulation. Temp. is measured by digital temp indicator.

The whole setup is mounted on a rigid M.S. frame.

Experimentation:-

To study the decomposition of Hydrogen Peroxide with water in presence of an Iodide catalyst.

To predict the degree of conversion from time-temperature data.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003