

Product Name :
Beakers

Product Code :
NLAB-EDUSCHOOLEQ122001



Description :

Beakers

Technical Specification :

Lab-Wares Beakers, moulded in Polypropylene, have excellent clarity and very good chemical resistance.

These autoclavable Beakers have prominently raised graduations on the outer surface.

Pouring a solutions becomes very convenient because of its tapered spout.

Material:- Polypropylene

Packing	
---------	--

50 Pts.	
100 Pts.	
250 Pts.	
500 Pts.	
600 mmml	
2500 ml	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003