

**Product Name :**  
Bed Side Locker-B

**Product Code :**  
NLABMEDICAL-M2730003

**Description :**

Bed Side Locker-B

**Technical Specification :**

The locker machine pressed from 20 gauge CRCA sheet

All three sides closed having one box and a drawer

Stainless Steel top

Mounted on 5 cm DIA. Swivelling castors

Finish in epoxy powder coated.

Bedside lockers provided with two storage cabinet and single draw under the top and space for keeping utilities

Bedside lockers optionally available without draw under the top

A high quality material is used for the Bedside Cabinets

These Bedside Lockers are helpful in keeping the medicines and other useful things in an organized manner

Moreover, these products are ideally placed at the bedside so that all the stored products remain within the reach

Side locker made with heavy duty CRCA MS sheet, duly machine pressed.

Provision of one small drawer

A big cabinet located below the Top

Stainless Steel top

Provision of one small drawer

A big cabinet located below the Top

Frame made of mild steel epoxy powder coated, top, front of the drawer

lower door made of laminated wood heat and stain resistance, mounted on four swivelling castors

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003