

<b>Product Name :</b> Bed Side Locker-B	<b>Product Code :</b> NLABMEDICAL-M2730003
<b>Description :</b>	
Bed Side Locker-B	
<b>Technical Specification :</b>	
<p>The locker machine pressed from 20 gauge CRCA sheet</p> <p>All three sides closed having one box and a drawer</p> <p>Stainless Steel top</p> <p>Mounted on 5 cm DIA. Swivelling castors</p> <p>Finish in epoxy powder coated.</p> <p>Bedside lockers provided with two storage cabinet and single draw under the top and space for keeping utilities</p> <p>Bedside lockers optionally available without draw under the top</p> <p>A high quality material is used for the Bedside Cabinets</p> <p>These Bedside Lockers are helpful in keeping the medicines and other useful things in an organized manner</p> <p>Moreover, these products are ideally placed at the bedside so that all the stored products remain within the reach</p> <p>Side locker made with heavy duty CRCA MS sheet, duly machine pressed.</p> <p>Provision of one small drawer</p> <p>A big cabinet located below the Top</p> <p>Stainless Steel top</p> <p>Provision of one small drawer</p> <p>A big cabinet located below the Top</p>	

Frame made of mild steel epoxy powder coated, top, front of the drawer

lower door made of laminated wood heat and stain resistance, mounted on four swivelling castors

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003