


<b>Product Name :</b> Bed Side Locker-C	<b>Product Code :</b> NLABMEDICAL-M2730004
--	---



**Description :**

Bed Side Locker-C

**Technical Specification :**

The locker box with deep drawer is made of commercial board flushed with laminate, assembled on square chrome plated tubes

The locker machine pressed from 20 gauge CRCA sheet

All three sides closed having one box and a drawer

Stainless Steel top

Mounted on 5 cm DIA. Swivelling castors

Finish in epoxy powder coated.

Bedside lockers provided with two storage cabinet and single draw under the top and space for keeping utilities

Bedside lockers optionally available without draw under the top

A high quality material is used for the Bedside Cabinets

These Bedside Lockers are helpful in keeping the medicines and other useful things in an organized manner

Moreover, these products are ideally placed at the bedside so that all the stored products remain within the reach

Side locker made with heavy duty CRCA MS sheet, duly machine pressed.

Size Top :(40 cm sq. x 82 H) cm

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003