

Product Name :
Bed Side Locker-D

Product Code :
NLABMEDICAL-M2730005

Description :

Bed Side Locker-D

Technical Specification :

Enamelled sheet steel frame

Full length body

One piece drawer

Laminated drawer face and top

Interior fitted with on steel shelf and perforated bottom

The locker box with deep drawer is made of commercial board flushed with laminate, assembled on square chrome plated tubes

The locker machine pressed from 20 gauge CRCA sheet

All three sides closed having one box and a drawer

Stainless Steel top

Mounted on 5 cm DIA. Swivelling castors

Finish in epoxy powder coated.

Bedside lockers provided with two storage cabinet and single draw under the top and space for keeping utilities

Size:- (40 cm sq. x 82 H) cm.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments