

Product Name :
Bedside Locker Deluxe

Product Code :
NLABMEDICAL-M2810001



Description :

Bedside Locker Deluxe

Technical Specification :

Frame Work of Bedside Locker Deluxe is made of mild steel pressed sheets around square tubes.

It has a stainless steel top with raised edges.

It is closed from three sides.

It comes with a cupboard box and one drawer for storage.

Bedside Locker Deluxe is Pre Treated Epoxy Powder Coated.

Product Size : 40 L x 40 W x 82 H cms.

We are Manufacturers, Exporters, and Suppliers of Bedside Locker Deluxe at tender and export competitive prices from India.

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

Naugralabequipments