

Product Name :
DIAC Characteristics Trainer Analog

Product Code :
NLAB-ELECTRICALAB210021



Description :

DIAC Characteristics Trainer Analog

Technical Specification :

Analog Meters:

Ammeter 2.5mA/25mA DC.

Voltmeter 50V DC.

Power supplies:

DC Supply IC Regulated 0-50V DC, 150mA.

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Components are mounted on the panels are:

DIAC DB3

Voltage Control through Potentiometer.

Salient Features:

Instruction manual.

Operated on Mains power 230V, 50Hz.

Patch Cords 4mm.

Supply provided through the 3 pin computer socket with appropriate fuse

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 10"x8"

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003