

**Product Name :**  
Digital Penetrometer

**Product Code :**  
NLAB-BTMENLAB150010



**Description :**

Digital Penetrometer

**Technical Specification :**

Centrifuge Extractor, Electrical Operation, Capacity 1500g, with a Dimmer stat for speed control from 2,400 to 3,600 rpm. Suitable for operation on 230 V, 50 Hz, Single Phase, A.C. supply.

Optional Extra: Filter Paper Discs, Set of 25

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003