

Product Name : Dynamic Trolley	Product Code : NLAB-EDUSCHOOLEQ1730215
Description :	
Dynamic Trolley	
Technical Specification :	
<p>This pair of identical wooden trolleys 300 mm long, mounted on three wheels.</p> <p>The high quality pair of identical 300mm long wooden trolleys have three low-friction wheels, a strong spring plunger.</p> <p>Supplied complete with 3 elastic rings, 3 springs, 4 pins (for stacking), 2 rubber corks, 2 needles and 1 releasing pin.</p> <p>It has a provision for stacking the trolleys on top of each other with removable pins.</p> <p>A spring loaded rod is fitted to give impulse when released.</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003