

**Product Name :**  
Educational Lab Copper Calorimeter

**Product Code :**  
NLAB-EDUSCHOOLEQ34051



**Description :**

Educational Lab Copper Calorimeter

**Technical Specification :**

A Copper Calorimeter for experimenting with the cooling of liquids.

The plastic lid has a built-in stirrer and a hole for a thermometer to record the temperature changes.

The calorimeter consists of a nickel-plated copper outer vessel with a lagged copper inner vessel.

Apparatus for experimenting with the cooling of liquids.

Nickel-plated copper outer vessel with a lagged copper inner vessel.

Outer vessel is 100mm high and 65mm in diameter.

Inner vessel is 75mm high and 50mm in diameter.

Built-in stirrer.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003