

Product Name :
Flow Through Helical Coil

Product Code :
NLAB-EDUCATIONTVET16018

Description :

Flow Through Helical Coil

Technical Specification :

Flow Through Helical Coil -

The whole set-up is well designed and arranged in a good quality painted structure

Helical Coil Material Stainless Steel Dia - 12"

Water Circulation FHP Pump, Crompton/Sharp make

Manometer 2 Nos

No. of Turns 15

Coil Tube Dia OD 21mm, ID 17mm (Approx.)

Sump Tank Material Stainless Steel/PVC, Capacity 50 Ltrs

Flow measurement Rotameter

Piping GI and PVC

Control panel comprising of Standard make on off switch, Mains Indicator etc

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003