

Product Name : Flow Through Helical Coil	Product Code : NLAB-EDUCATIONTVET16018
Description : Flow Through Helical Coil	
Technical Specification : Flow Through Helical Coil - The whole set-up is well designed and arranged in a good quality painted structure Helical Coil Material Stainless Steel Dia - 12" Water Circulation FHP Pump, Crompton/Sharp make Manometer 2 Nos No. of Turns 15 Coil Tube Dia OD 21mm, ID 17mm (Approx.)	

Sump Tank Material Stainless Steel/PVC, Capacity 50 Ltrs

Flow measurement Rotameter

Piping GI and PVC

Control panel comprising of Standard make on off switch, Mains Indicator etc

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003