

Product Name :
Fowler Bed Super

Product Code :
NLABMEDICAL-M2800002

Description :

Fowler Bed Super India, Fowler Bed Super manufacturer India, suppliers, Medical Fowler Bed Super India, Hospital Fowler Bed Super India

Technical Specification :

Backrest adjustable on crank mechanism

Stainless steel handle for backrest adjustment

These fowler beds are equipped with easy lifting or collapsible side rails to ensure patient safety

Swivelling castors with two brakes

Stainless Steel or ABS is used for the head and foot panels

Four locations for saline rod & safety side railing Finish in epoxy powder coated

Four section top of uniformly perforated CRCA sheet

Backrest and knee rest section maneuvered by two crank mechanisms from the foot end

ABS bows easily removable for access to the patient

Trolley type base mounted on 125 mm die

Frame work of precise steel tubes

Size: (198 L x 90 W x 60 H) cm.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

