

**Product Name :**  
Fowler Bed

**Product Code :**  
NLABMEDICAL-M2790005

**Description :**

Fowler Bed India, Fowler Bed manufacturer India, suppliers, Medical Fowler Bed India, Hospital Fowler Bed India

**Technical Specification :**

To ensure its quality, this bed is strictly tested on various defined quality parameters by our quality controllers.

Fowler Beds made of high quality materials, components and accessories

Fowler Beds with adjustable back section and knee-rest

High quality castors

Backrest adjustable on crank mechanism

Stainless steel handle for backrest adjustment

These fowler beds are equipped with easy lifting or collapsible side rails to ensure patient safety

Stainless Steel or ABS is used for the head and foot panels

Four section top of uniformly perforated CRCA sheet.

Backrest and knee rest section maneuvered by two crank mechanisms from the foot end.

ABS bows easily removable for access to the patient

Trolley type base mounted on 125 mm die

General Specification:-

Overall Size: Mattress size

{LxWxH}: {LxW}

{2160x1015x600}: {1980x910}

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003