

Product Name : Fowler Bed	Product Code : NLABMEDICAL-M2790005
Description : Fowler Bed India, Fowler Bed manufacturer India, suppliers, Medical Fowler Bed India, Hospital Fowler Bed India	
Technical Specification : To ensure its quality, this bed is strictly tested on various defined quality parameters by our quality controllers. Fowler Beds made of high quality materials, components and accessories Fowler Beds with adjustable back section and knee-rest High quality castors Backrest adjustable on crank mechanism Stainless steel handle for backrest adjustment These fowler beds are equipped with easy lifting or collapsible side rails to ensure patient safety Stainless Steel or ABS is used for the head and foot panels Four section top of uniformly perforated CRCA sheet. Backrest and knee rest section maneuvered by two crank mechanisms from the foot end. ABS bows easily removable for access to the patient Trolley type base mounted on 125 mm die General Specification:- Overall Size: Mattress size {LxWxH}: {LxW} {2160x1015x600}: {1980x910}	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003