

**Product Name :**  
Hydraulic Valves and Fittings Test Stand

**Product Code :**  
TK198



**Description :**

Hydraulic Valves and Fittings Test Stand

**Technical Specification :**

Test stand on which to mount industrial valves and fittings  
Pressure testing of valves and fittings  
Hand-operated piston pump to generate the test pressure, a return valve to relieve the system pressure, and a manometer for pressure measurement  
2 different sizes of mounting flange with blank flange and flange seal  
Connection of pump and test flange via pressure hose  
Test medium: water  
Mobile frame with collector tray and ball valve to drain  
Water storage tank  
Part of the Practice Line for assembly, maintenance and repair

**Technical Data:**

Piston pump with tank  
test pressure: 0...60bar  
tank capacity: 12L  
manometer: 0...60bar  
Mounting flanges for valves and fittings under test

**Dimensions and weight:**

Length x Width x Height: 1000x750x1200mm

Weight: 80kg

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003