

<b>Product Name :</b> ICU Bed Manual	<b>Product Code :</b> NLABMEDICAL-M2790008
---	---

--

**Description :**

ICU Bed Manual India, ICU Bed Manual manufacturer India, suppliers, Medical ICU Bed Manual India, Hospital ICU Bed Manual India

**Technical Specification :**

These beds are made in strict compliance with the international standards

Top of 1.25mm thick CRCA sheet perforated uniformly

Base and mainframe work of precise steel tubes

Base mounted on 125cm DIA. swiveling castors, two with break.

The backrest and knee rest as also trendelenberg and reverse trendelenberg obtainable on smooth crank mechanism

The height of the Bed can be adjusted to facilitate close observation by the Medical staff

Stainless steel bows with laminated panels, easily removable for access to the patient

Fitting: SS saline rod with four locations, protective SS railing. Epoxy powder coated finish

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003