

**Product Name :**  
Instrumentation Amplifier Using OP-AMP Trainer

**Product Code :**  
NLAB-ELECTRONICSAB110018



**Description :**

Instrumentation Amplifier Using OP-AMP Trainer

**Technical Specification :**

**SCOPE OF LEARNING:**

Study of Instrumentation Amplifier

**TECHNICAL SPECIFICATIONS:**

**Power Supplies:**

DC Power Supply IC Regulated 0-5V DC, 150mA. (2Nos.)

Operated on Mains power 230V, 50Hz +10%

**Digital Meters:**

Voltmeter 20V

Components are mounted on the panels are:

Resistors

IC LM741

Variable Resistor For Gain Setting

SALIENT FEATURES:

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Instruction manual.

Connections are brought out through 2mm Colored Sockets.

Patch Cords 2mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003