

Product Name :
Kinetic Theory Model

Product Code :
NLAB-EDUSCHOOLEQ174055



Description :

Kinetic Theory Model

Technical Specification :

To simulate the behaviour of a gas as its molecules gain energy. Comprises a 50 mm diam. clear acrylic tube of 380 mm height, its lower end enclosing a vibrator beneath a diaphragm. With a loose cover for the upper end, two expanded polystyrene pistons and a supply of metal spheres. Complete with a 6V electric motor to drive the vibrator.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003