

Product Name : Kinetics of Dissolution of Benzoic Acid	Product Code : NLAB-EDUCATIONTVET16033
Description : Kinetics of Dissolution of Benzoic Acid	
Technical Specification : Kinetics of Dissolution of Benzoic Acid - The setup consists of a Reactor in which solid cylinder of Benzoic Acid is placed and water is stirred by means of a variable speed mixer. Baffles are also provided. A Peltizer is provided to make the solid cylinder size 1? diameter and 2? height of Benzoic Acid. The Reactor is made of Stainless Steel. Water Supply and Drain. Electricity Supply: 1 Phase, 220 V AC, 0.5 kW. To study the Kinetics of dissolution of Benzoic Acid.	

Instruments, Laboratory Glassware, and Chemicals required for analysis.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003