

**Product Name :**  
Kinetics of Dissolution of Benzoic Acid

**Product Code :**  
NLAB-EDUCATIONTVET16033

**Description :**

Kinetics of Dissolution of Benzoic Acid

**Technical Specification :**

Kinetics of Dissolution of Benzoic Acid -

The setup consists of a Reactor in which solid cylinder of Benzoic Acid is placed and water is stirred by means of a variable speed mixer. Baffles are also provided.

A Peltizer is provided to make the solid cylinder size 1? diameter and 2? height of Benzoic Acid.

The Reactor is made of Stainless Steel.

Water Supply and Drain.

Electricity Supply: 1 Phase, 220 V AC, 0.5 kW.

To study the Kinetics of dissolution of Benzoic Acid.

Instruments, Laboratory Glassware, and Chemicals required for analysis.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003