

**Product Name :**  
Lockable Box for Storage

**Product Code :**  
NLAB-SCHOOLABS25013



**Description :**

Lockable Box for Storage

**Technical Specification :**

Metal storage box, 790x590x245mm, lockable, stackable, supplied with two code padlocks

Technical Specifications:

Metal Box:

Made of high density, durable lightweight aluminum. Robust, lockable, stackable. Water proof, dustproof and impact resistant.

Material, recycled aluminum, thickness 1mm.

- External dimensions: 790x590x245mm
- Internal dimensions: 764x574x220mm
- 2 Locking front side lock for code padlocks
- 2 Strap for cover nylon

Code Padlocks, set of 2, 2 sets are provided with the box:

The padlock consists of shackle, a dial, locking parts, and means for combination changing and combination scrambling when locked.

Locks are fixed inside the cover of the box.

Weight and Volume:

Volume: 0.12cbm

Weight: 6.5kg/Weight pr. gsm: 135kg

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003