

Product Name :
Medicine Cabinet

Product Code :
NLABMEDICAL-M2730008

Description :

Medicine Cabinet

Technical Specification :

Medicine Cabinets are usually used in laboratories, hospital, research centers for keeping the drugs, medicines and lab chemicals safe

Made of Mild Steel duly epoxy powder coated

Top section with glass doors and two shelves

Middle section with four drawers

Lower section with two doors and one shelve with lock.

These are manufactured in various designs such as with glass doors or without glass doors, with standard specifications or as per customer's specifications

Frame Work Constructed with Heavy duty Ms Steel Sheet

Two Glass Doors Available With Individual Key Locking

Three Adjustable Glass Selves

One Bottom Cabinet with Key Locking Facility

Epoxy coated mild steel framework.

Stainless steel bottom

Size :-Top: {25 x 90 x 85 cm (H)}, {Base: 45 x 90 x 92 cm (H)}

Specifications:

Design: Frame

Feature: Eco-Friendly, Rust Proof, Hygienic

General Use: Indoor Furniture

Furniture Type: Laboratory Furniture

Laboratory Furniture Type: Laboratory Cabinets

Material: Stain steel

Metals Type: Steel

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003