

Product Name : Medicine Cabinet	Product Code : NLABMEDICAL-M2730008
Description :	
Medicine Cabinet	
Technical Specification :	
<p>Medicine Cabinets are usually used in laboratories, hospital, research centers for keeping the drugs, medicines and lab chemicals safe</p> <p>Made of Mild Steel duly epoxy powder coated</p> <p>Top section with glass doors and two shelves</p> <p>Middle section with four drawers</p> <p>Lower section with two doors and one shelf with lock.</p> <p>These are manufactured in various designs such as with glass doors or without glass doors, with standard specifications or as per customer's specifications</p> <p>Frame Work Constructed with Heavy duty Ms Steel Sheet</p> <p>Two Glass Doors Available With Individual Key Locking</p> <p>Three Adjustable Glass Selves</p> <p>One Bottom Cabinet with Key Locking Facility</p> <p>Epoxy coated mild steel framework.</p> <p>Stainless steel bottom</p> <p>Size :-Top: {25 x 90 x 85 cm (H)}, {Base: 45 x 90 x 92 cm (H)}</p> <p>Specifications:</p> <p>Design: Frame</p>	

Feature: Eco-Friendly, Rust Proof, Hygienic

General Use: Indoor Furniture

Furniture Type: Laboratory Furniture

Laboratory Furniture Type: Laboratory Cabinets

Material: Stain steel

Metals Type: Steel

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003