

Product Name :
Mosfet Characteristics Trainer Analog

Product Code :
NLAB-ELECTRICALLAB210037



Description :

Mosfet Characteristics Trainer Analog

Technical Specification :

Analog Meters:

Voltmeter 5V DC.

Ammeter 30mA DC.

Voltmeter 30V DC.

Power supplies:

DC Supply IC Regulated 0-5V DC, 150mA.

DC Supply IC Regulated 0-30V DC, 150mA.

Front panel built with high class insulated Printed Circuit Board sheet with well printed circuits and symbols.

Fuse for Short Circuit protection

Components are mounted on the panels are:

Mosfet

Voltage Control through Potentiometer.

Salient Features:

Instruction manual.

Operated on Mains power 230V, 50Hz.

Patch Cords 4mm.

Supply provided through the 3 pin computer socket with appropriate fuse

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x8"

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003