

Product Name :
Multi Therm Cool Heat Shaker

Product Code :
NLAB-EDUSCHOOLEQ45018

Description :

Multi Therm Cool Heat Shaker

Technical Specification :

The Multi Therm shaker is a temperature controlled vortexer, useful for a variety of molecular biology applications.

Integral over-temperature control ensures long life, safety, and sample integrity.

Speed, time and temperature settings are continuously visible on the LCD, simultaneously showing both actual and selected values.

Despite its modest 8 x 11 in. footprint, all popular tubes, and plates are compatible, including as many as six 50ml tubes.

Speed, 200 to 1500 rpm (or off)

Orbit:

Temp. Range: 20 to 100°C (H5000-HC)
Ambient +5 to 100°C (H5000-H)

Temp. Increment:

Temp. Accuracy:

Cap. Selected Block

Speed Increment:

5.000 g per Rate in.

2.500 g per Rate in.

Dimensions (W x D x H) 21 x 29 x 21 cm

Weight (kg)

Timer: to 99 hours 99min. / continuous

Electrical 230V, 50-60Hz, 500W

Warranty:

Speed, 200 to 1500 rpm (or off)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003