

Product Name : My Bath Series Water Bath	Product Code : NLAB-EDUSCHOOLEQ45020
Description :	
My Bath Series Water Bath	
Technical Specification :	
<p>The myBath series is designed to guarantee superior accuracy and performance in a compact footprint.</p> <p>Adjustable in 0.1°C increments, the temperature is shown on a large, 3-digit display.</p> <p>High efficiency insulation maintains temperature stability and uniformity within the bath, while keeping the housing cool to the touch.</p> <p>With four sizes to choose from, there is a myBath to fit every lab.</p> <p>The 4, 8 and 12L models include easy to use controls allow the user to digitaly select and monitor settings.</p> <p>All myBath models feature seamless, stainless steel chambers for easy cleaning and long life.</p> <p>Racks are available to accommodate 1.5/2.0ml, 15ml and 50 ml tubes within the 4L, 8L and 12L baths.</p> <p>The myBath 4L is supplied with all three racks.</p> <p>The chambers resist rust and have no corners to collect dirt or contamination.</p> <p>The larger 8L and 12L units have drains for hassle-free water changes.</p> <p>The 4L, 8L, and 12L models include a see-through cover that can be removed, or flipped back on a rear hinge to conserve bench space.</p> <p>The lid's gabled design will divert condensation to the front and back of the bath, away from samples.</p> <p>Specification:-</p>	

Quik-CAL user calibration function.

Mybath 4L Includes racks for 1.5/2.0, 15 and 50ml tubes.

Seamless, stainless steel chamber.

Hinged lid (included) also removable.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003