

Product Name :
My Bath Series Water Bath

Product Code :
NLAB-EDUSCHOOLEQ45020

Description :

My Bath Series Water Bath

Technical Specification :

The myBath series is designed to guarantee superior accuracy and performance in a compact footprint.

Adjustable in 0.1°C increments, the temperature is shown on a large, 3-digit display.

High efficiency insulation maintains temperature stability and uniformity within the bath, while keeping the housing cool to the touch.

With four sizes to choose from, there is a myBath to fit every lab.

The 4, 8 and 12L models include easy to use controls allow the user to digitaly select and monitor settings.

All myBath models feature seamless, stainless steel chambers for easy cleaning and long life.

Racks are available to accommodate 1.5/2.0ml, 15ml and 50 ml tubes within the 4L, 8L and 12L baths.

The myBath 4L is supplied with all three racks.

The chambers resist rust and have no corners to collect dirt or contamination.

The larger 8L and 12L units have drains for hassle-free water changes.

The 4L, 8L, and 12L models include a see-through cover that can be removed, or flipped back on a rear hinge to conserve bench space.

The lid's gabled design will divert condensation to the front and back of the bath, away from samples.

Specification:-

Quik-CAL user calibration function.

Mybath 4L Includes racks for 1.5/2.0, 15 and 50ml tubes.

Seamless, stainless steel chamber.

Hinged lid (included) also removable.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003