

Product Name :
Oven For Kiln Drying

Product Code :
NLAB-MACHINELABWK33010



Description :

Oven For Kiln Drying

Technical Specification :

Made up of fixed frame capable of holding about 500 m3 of sawn timber

Rolling superstructure

With circulating fans of airspeed 4.5m/s

Kiln temperature sensors up to 80oC

Come complete with supplementary heating option of hybrid gas heating system

Type Solar Kiln

Wood Workshop Portable Tools Oven For Kiln Drying Manufacturer,Supplier and Oven For Kiln Drying Exporter in India

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003