

Product Name :
Overbed Bedside Table Adjustable Gear Handle

Product Code :
NLABMEDICAL-M2690001



Description :

Overbed Bedside Table Adjustable Gear Handle

Technical Specification :

The frame is made of Rectangle & Square Tube.

The body section is mounted on 5cms castors.

The size of the Laminated top is 75L x 40W cms.

The height of the Overbed Bedside table is adjustable manually by Gear Handle.

Knock Down Construction to save freight.

We are Manufacturer, Exporters, and Suppliers of Overbed Bedside Table Adjustable Gear Handle at tender and export competitive prices from India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

