

**Product Name :**  
Overbed Bedside Table Adjustable Manual

**Product Code :**  
NLABMEDICAL-M2690003



**Description :**

Overbed Bedside Table Adjustable Manual

**Technical Specification :**

Specification :-

Frame is made of Rectangle & Square Tube.

The body section is mounted on 5cms castors.

The size of Laminated top is 75L x 40W cms.

The height of Overbed Bedside table is adjustable manually by side Knob.

Knock Down Construction to save freight.

We are Manufacturer, Exporters and Suppliers of Overbed Bedside Table Adjustable Manual at tender and export competitive prices from India.

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## Naugralabequipments