

Product Name :
Plasma Extractor

Product Code :
NLAB-EDUSCHOOLEQ46003



Description :

Plasma Extractor

Technical Specification :

The Plasma Extractor is designed to express.

Blood components from blood collection bags.

The Main bag is to be placed between the fixed stainless steel and movable transparent plate where the power spring load provides pressure to transfer the fluid into designated satellite bag.

External size (W*D*H): 165*250*250 (mm)

Gross weight 2.8kg

Clamping Force 20N

Package size (W*D*H): 310*310*240 (mm)

Applicable Blood Bag 100ml, 200ml, 300ml, 400ml

Net weight 2.5kg

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003