

Product Name :
Reverse Osmosis and Ultra Filtration Pilot Plant India

Product Code :
NLAB-ENGINEERINGLB36036



Description :

Reverse Osmosis and Ultra Filtration Pilot Plant

Technical Specification :

Safety pressure switches.

Heat exchanger

4x thermocouples with PT100 probes protected by stainless steel

shield sleeves

2x bourdon gauges

2x variable area flow meters for permeate and concentrate.

Safety valve

Electrical control unit

Emergency stop button.

Feed tank with a capacity of 30 liters.

Tank for collecting the permeate with a capacity of 30 liters.

Piston pump.

2x membranes with typical permeate flow rate of 5 to 50 milli liters per minute

Power Supply:-

400 Vac 50 Hz three-phase -3000 VA

We are well-known manufacturers, OEM suppliers of Reverse Osmosis and Ultra Filtration Pilot Plant for Chemical Engineering Lab Equipments. Contact us for high quality Reverse Osmosis and Ultra Filtration Pilot Plant for Chemical Engineering Lab Equipments for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003