

Product Name : Saline Stand	Product Code : NLABMEDICAL-M2730010
Description :	
Saline Stand	
Technical Specification :	
Strong Stainless Steel construction mounted on Moulded Base fitted with Five swivel castors, 50mm dia. Chrome plated rod with double hooks height adjustable from approx 1220 to 2340 mm. Also Available With Stainless Steel Base Features: Heavy duty castor wheels Rust proof baseLightweight Highly durable Color White Product Details: Finished: Epoxy, powder coated Material: Stainless steel Base: Heavy Duty Plastic	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

