


<b>Product Name :</b> Saline Stand	<b>Product Code :</b> NLABMEDICAL-M2730010
---------------------------------------	---



**Description :**

Saline Stand

**Technical Specification :**

Strong Stainless Steel construction mounted on Moulded Base fitted with Five swivel castors, 50mm dia.

Chrome plated rod with double hooks height adjustable from approx 1220 to 2340 mm.

Also Available With Stainless Steel Base

Features:

Heavy duty castor wheels

Rust proof baseLightweight

Highly durable

Color White

Product Details:

Finished: Epoxy, powder coated

Material: Stainless steel

Base: Heavy Duty Plastic

**Naugralabequipments**  
**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003

