

Product Name :
Surgeons Stool

Product Code :
NLABMEDICAL-M2730002



Description :

Surgeons Stool

Technical Specification :

Mounted on five molded five swiveling castors, strong steel construction, top and back upholstered with cushioned washable material.

Height of the stool is adjustable by Gas Pump controlled by hand lever at lower side of the seat.

Size: 300mm sq x H 510 mm

Features:

Adjustable height

Cushioned top

Superior finish

Description:

CNC machine pressed & bent five prong MS frame work mounted on 50mm dia. non-rusting castor wheels

Foot ring is of approximate dia 485 mm

Finish:

Mechanically pre-treated and epoxy powder coated.

Optional:

Foot ring of Stainless steel

Specifications:

Overall Approximate size: 385 mm Dia

Height is adjustable by imported gas spring from 500 mm to 650 mm

Cushioned top of approximate 385 mm dia

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003