

**Product Name :**  
TDS Meter Bench Top

**Product Code :**  
NLAB-EDUSCHOOLEQ46008



**Description :**

TDS Meter Bench Top

**Technical Specification :**

We own a broad range of reliable and accurate TDS Meters.

The main function of TDS Meters is to measure the Conductivity/TDS/Salinity of aqueous solutions.

Technical Specifications:-

Display: 3½ Digit LCD

Range: 2 Ranges 0-2000 ppm 0-20 ppt

Temperature Compensation: NA

Cell Constraint Adjustment: Pre-adjustable

Conductivity Cell: Platinum DIP

Measuring Frequency: 1000 Hz

Accuracy: ±1% FSI

Power: 1 x 9V Battery

Dimensions: 85 x 28 x 165 mm

Weight: 0.8 kg (approx.)

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003