

Product Name :
TDS Meter Bench Top

Product Code :
NLAB-EDUSCHOOLEQ46008



Description :

TDS Meter Bench Top

Technical Specification :

We own a broad range of reliable and accurate TDS Meters.

The main function of TDS Meters is to measure the Conductivity/TDS/Salinity of aqueous solutions.

Technical Specifications:-

Display: 3½ Digit LCD

Range: 2 Ranges 0-2000 ppm 0-20 ppt

Temperature Compensation: NA

Cell Constraint Adjustment: Pre-adjustable

Conductivity Cell: Platinum DIP

Measuring Frequency: 1000 Hz

Accuracy: $\pm 1\%$ FSI

Power: 1 x 9V Battery

Dimensions: 85 x 28 x 165 mm

Weight: 0.8 kg (approx.)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003