

**Product Name :**  
Triangular Chain Pulley Block

**Product Code :**  
NLAB-MECHANICALAB240052



**Description :**

Triangular Chain Pulley Block

**Technical Specification :**

We have been able to offer a wide range of superior quality Triangular Chain Pulley Blocks. These blocks are widely used for different medium duty applications for lifting purpose. The offered blocks are also tested on various parameters so as to offer a defect free range to our clients.

Featured:

Durability

Flexibility

High performance

Robust construction

Easy maintenance.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003